



**THT Biomaterials GmbH**  
extracellular platform technology  
**The Human Touch**

Office  
Schiffmühlenstraße 106/16  
1220 Vienna  
Fax +43 1 890 64 61  
Mobile: +43 676 93 66505  
Lab  
Vienna Biocenter (VBC6)  
Dr.-Bohr-Gasse 7  
1030 Vienna  
tht-biomaterials.com  
office@tht-biomaterials.com

## Product information

### **HUMAN PLACENTA Collagen-I (COL1), liquid**

Comprising highly purified human atelocollagen type 1 in 10 mM HCl, prepared from placenta tissue of individuals that have been shown by certified tests to be negative for antibodies to HIV, HEP-B and THPA (syphilis).

**Catalog number** THT0101

### **Product description**

Pepsin-solubilized atelocollagen allows growth and cultivation of many cell types. COL1 is derived from human origin, liquid and prepared at a concentration of 2 mg/mL (> 95% purity by SDS-PAGE).

### **Precautions and Disclaimer**

This product is only for R&D use. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling procedures.

### **Storage/Stability**

Delivery in cool packs. Following reconstitution, aliquot and freeze. Store COL1 stock solutions at -20 °C for up to 12 months.

### **Procedure**

*Note: The optimal concentration for cell attachment and culture may differ for different cell types, and experimentation may be required to determine the optimal conditions for your personal cell culture settings.*

### **Guideline for 2D coating tissue culture plates**

1. Add sufficient volume of COL1 to provide your desired coating concentration. If necessary, dilute your COL1 stock with 2 mM HCl.

A coating concentration of at least 0.1 µg/cm<sup>2</sup> is recommended, depending on the cell type used and your experimental settings.

It is important that the volume added to the well is sufficient to cover the whole growth surface.

2. Keep the plate completely covered and incubate for 60 min at room temperature.
3. Tilt the plate to allow excess COL1 to drain to the lowest point and remove the remaining excess material with a sterile pipette.
4. Air dry the plate, and use.

### **Guidelines for 3D gelling**

1. Dilute COL1 to your desired gel concentration using 2 mM HCl.
2. Mix 9 parts COL1 and 1 part 10X buffer or 10X cell culture media.
3. Add the mixture to your tissue culture plate.
4. Incubate the plate for 30 to 60 min at 37°C, and use. Gels must be handled carefully.

### **References**

1. Hackethal J, Mühleder S, Hofer A, Schneider KH, Prüller J, Hennerbichler S, Redl H, Teuschl A. *An Effective Method of Atelocollagen Type 1/3 Isolation from Human Placenta and Its In Vitro Characterization in Two-Dimensional and Three-Dimensional Cell Culture Applications*. Tissue Eng Part C Methods. 2017 May;23(5):274-285
2. Hackethal J, Hofer A, Hennerbichler S, Redl H, Teuschl A. *A comparison of enzymatic and non-enzymatic strategies to isolate extracellular matrix (ECM) proteins from human placenta and liposuction fat*. ALTEX Proceedings 8(1), p65, 2019.