


Technical Data Sheet

Part Number: LQ-632111

Description: LabQ Thin Wall 8- Strip PCR Tubes with individually

optically Flat caps, for QPCR flat optically clear cap,

Specifications



www.labconsulting.at Phone: +43 1244 28 55-56 office@labconsulting.at



Technical Data Sheet

Storage and Usability

Storage	Protect from direct sunlight and UV light
	Store dry at room temperature.
Operational temperature	-20°C to +105°C
Best by date	Five (5) years from date of manufacture
Usage	Designed for single use in PCR applications.
	Intended for research purposes only.

Endotoxin (Pyrogen)Testing: The product is tested to be in accordance with ANSI/AAMI ST 72:2002 | (R)2010 *Bacterial endotoxins test, Methodologies, routine monitoring and alternatives to batch testing.*

RNase/DNase Testing: This product is tested and is free of any detectable RNase/DNase contamination.

Bovine Spongiform Encephalopathy and Transmissible Spongiform Encephalopathy: This product complies with the latest revision of EMEA/401/01 "Note for Guidance on minimizing the risk of transmitting animal spongiform encephalopathy agents via human veterinary medicine products" by virtue of all bovine derived material having been processed per specific conditions of section 6.4 of EMEA/401/01

Composition

Color	Clear, have optically clear caps for qPCR Medical-grade virgin
-------	--

Phone: +43 1244 28 55-56

Material polypropylene